




Fast Respiration Meter Comparison Matrix	FRM S	FRM L	FRM XL
			
Used by*	<i>Small packhouse</i>	<i>Small and/or Large packhouse</i>	<i>Small&amp;Large packhouse / Film Converter</i>
Used for	<i>Cut vegetables, cut fruit, small whole fruit</i>	<i>Cut vegetables, cut fruit, small whole fruit</i>	<i>Cut vegetables, whole vegetable head, cut fruit, all size whole fruit, flowers</i>
Size Vessel	<i>Medium (7 liter)</i>	<i>Medium (7 liter)</i>	<i>Large (30 liter)</i>
Gas Flush possible	<i>NO</i>	<i>YES</i>	<i>YES</i>
Internal Temperature Control	<i>NO**</i>	<i>YES</i>	<i>YES</i>
Internet Connection	<i>NO</i>	<i>NO</i>	<i>YES</i>
Automatic upload of data	<i>NO</i>	<i>NO</i>	<i>YES</i>
Manual upload of data	<i>YES</i>	<i>YES</i>	<i>YES</i>
Remote Support	<i>NO</i>	<i>NO</i>	<i>YES</i>
Dimensions (LxWxH) [mm]	<i>602x668x460</i>	<i>620x635x1330</i>	<i>682x1126x1330</i>
Weight [kg]	<i>47</i>	<i>100</i>	<i>150</i>
To be placed on	<i>Table</i>	<i>Floor (includes wheels)</i>	<i>Floor (includes wheels)</i>
Power Consumption Max. [W]	<i>350</i>	<i>1900</i>	<i>1900</i>
Voltage [V]	<i>115/230V</i>	<i>115/230V</i>	<i>115/230V</i>
Release date	<i>February 2014</i>	<i>February 2010</i>	<i>February 2012</i>
<p>* PerfoTec advice based on common use</p> <p>**FRM needs to be placed in room with controlled temperature in order to properly measure respiration rates</p>			